

# LGPv 1420

## MediLine



|                                      |                                    |
|--------------------------------------|------------------------------------|
| <b>Total gross capacity</b>          | 1427 litre                         |
| <b>Ext. dimensions</b> in mm (w/d/h) | 1430/830/2150                      |
| <b>Int. dimensions</b> in mm (w/d/h) | 1250/700/1550                      |
| <b>Energy efficiency class</b>       |                                    |
| Energieverbrauch in 365 Tagen        | 3244 kWh                           |
| Climate rating <sup>1</sup>          | 5                                  |
| Cooling system                       | –                                  |
| Defrosting                           | automatic                          |
| Temperature range                    | -10°C to -26°C                     |
| Gastronorm                           | - mm                               |
| Housing material/colour              |                                    |
| Door/lid material                    | Steel                              |
| Liner Material                       | 304-grade stainless steel          |
| Type of control                      | electronic with text display       |
| Temperature display                  | external digital                   |
| Interior light                       | –                                  |
| Fault: warning signal                | visual and audible                 |
| Interface                            | RS 485                             |
| volt-free alarm contact              | yes                                |
| Shelves                              | 8                                  |
| Adjustable Shelving                  | –                                  |
| Number of drawers                    | –                                  |
| Number of baskets                    | –                                  |
| Maximum number of baskets            | –                                  |
| Shelf material                       | grid shelves plastic-coated        |
| Shelf loading                        | 60 kg                              |
| Height of adjustable legs            | -                                  |
| Castors/type of castors              | yes/swivel castors                 |
| Door handle                          | full-length integrated door handle |
| Foot pedal opener                    | –                                  |
| Self-closing door                    | yes                                |
| Insulation                           | 83 - 60 mm                         |
| Defrost water drain                  | –                                  |
| Refrigerant                          | R 290                              |
| Airborne acoustical noise emissions  | 60 dB(A)                           |
| Gross/net weight                     | 262,0/219,0 kg                     |
| EAN-Nr.                              | 9005382137215                      |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |
|                                      |                                    |

Subject to modification. Timestamp: 03.11.2016

# LIEBHERR



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.