

# LGUex 1500

## MediLine



<b>Total gross capacity</b>	139 litre
<b>Ext. dimensions</b> in mm (w/d/h)	600/615/820
<b>Int. dimensions</b> in mm (w/d/h)	454/450/663
<b>Energy consumption</b> in 24 h	0,926 kWh
Climate rating <sup>1</sup>	SN-ST
Cooling system	static
Defrosting	manual
Temperature range	-9°C to -26°C
Gastronorm	-
Housing material/colour	steel/white
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	electronic controller
Temperature display	external digital
Interior light	-
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	4
Adjustable Shelving	-
Number of drawers	3
Number of baskets	4
Maximum number of baskets	-
Shelf material	evaporator plates
Shelf loading	24 kg
Height of adjustable legs	-
Castors/type of castors	-
Door handle	Ergonomic door handle
Foot pedal opener	-
Self-closing door	yes
Insulation	55 - 55 mm
Defrost water drain	-
Refrigerant	R 600a
Airborne acoustical noise emissions	45 dB(A)
Gross/net weight	44/42 kg
EAN-Nr.	9005382171011

Subject to modification. Timestamp: 22.03.2016

# LIEBHERR



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

