


Life Science Tools & Technology

UV Products | Light Sources | Lab Products



UV Blocking Eyewear



Eyewear protects the eyes from harmful UV. Goggles and Faceshield comply with ANSI-Z87.1 1986 specifications for eye and face protection devices.

Eye and face protection is essential for anyone working with ultraviolet sources which can cause burning to unprotected eyes and skin. Eyewear blocks wavelengths from 200 - 400 nm.

- **Spectacles** are made of impact resistant polycarbonate. Top and side-shield design provides UV protection for the eyes. Orange spectacles block blue light that can cause haze
- **Goggles** are constructed of specially-formulated plastic for optimum viewing contrast and lessened eye fatigue. Goggles feature a lightweight design and flexible rubber trim. Air vents allow air flow while blocking ultraviolet transmission
- **Faceshield** covers the face and neck area

Order Number Description

	UV Blocking Eyewear
98-0002-01	UVC-303 Spectacles; 2H x 6W in. (51 x 152mm)
38-0240-01	UVC-310 Spectacles, orange; 2H x 6W in. (51 x 152mm)
98-0002-02	Goggles UVC-503; 3H x 6W in. (76 x 152mm)
98-0002-03	Replacement Goggle Lens
98-0002-04	Faceshield UVC-803; 17dia x 8.25H in. (432 x 210mm)
38-0151-01	Replacement Shield



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

