



Argos RotoFlex™

Benchtop Tube Rotator designed for a variety of applications!

The RotoFlex will accommodate a variety of tubes from 0.5 - 50mL by simply changing the rotisserie. No tools required! The mixing angle can be adjusted from 0-90° and each of the 4 tube panels can be aligned independently. The 360° rotation ensures thorough mixing of all samples at the same time while rotating at 22 rpm.

To change the rotation direction - simply touch the rotisserie paddle while in motion and the RotoFlex automatically begins to rotate in the opposite direction.

RotoFlex can be used in incubators & cold room settings.
CE Approved / RoHS Compliant.

Cat. No.	Description	Qty.
R2000	RotoFlex 120V	1/EA
R2001	RotoFlex 220V	1/EA
R2002	Rotisserie for .5 - .8mL tubes	2/PK
R2003	Rotisserie for 1.5 - 2mL tubes	2/PK
R2004	Rotisserie for 5 - 7mL tubes	2/PK
R2005	Rotisserie for 10 - 15mL tubes	2/PK
R2006	Rotisserie for 50mL tubes	1/PK

**Unit is supplied with base unit and rotisseries
(1 pack of each size).**

Ideal for:

- Hybridization
- DNA Purification
- Sample mixing / resuspension
- Protein Binding
- Immunoprecipitation



Removable Inserts - easily snap in and out of place



Change inserts to accommodate varying tube sizes



Simply touch to engage auto reverse feature



New Item!

ARGOS

TECHNOLOGIES

| Distributed by: |



LIT. NO. 9950L-A / REV0112



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

