



Function / construction:

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **easy alignment:** adjusting aids to compensate for uneven floor
- **ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

Available equipment

- 2 shelves (height 1105 mm) or 3 shelves (height 1950 mm) galvanised with bottom collecting sump
- doors with glass cutout – the hazardous materials can be seen at a glance, even when the doors are closed

Cabinet for chemicals

CS-CLASSIC-G model CS.195.105.WDFW:0003

Body colour Anthracite grey (RAL 7016), with wing doors and glass cutout in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30529



top
price
**AVAILABLE
FROM STOCK
IMMEDIATELY!**




Cabinet for chemicals

CS-CLASSIC-G model CS.110.105.WDFW:0002

Body colour Anthracite grey (RAL 7016), with wing doors and glass cutout in Traffic red (RAL 3020), interior equipment with 2 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30527



Function / construction:

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **easy alignment:** adjusting aids to compensate for uneven floor
- **ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

Available equipment

- 2 shelves (height 1105 mm) or 3 shelves (height 1950 mm) galvanised with bottom collecting sump

Cabinet for chemicals

CS-CLASSIC model CS.195.105:0003

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30528



top
price

**AVAILABLE
FROM STOCK
IMMEDIATELY!**



Cabinet for chemicals

CS-CLASSIC model CS.110.105:0002

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 2 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30526

C-LINE



Safe and proper storage of non flammable chemicals and toxic substances in work spaces

Function / construction:

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **easy alignment:** adjusting aids to compensate for uneven floor

- **ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

Available equipment

- shelves – optionally in galvanised or painted version, height adjustable
- perforated insert useable as storage area
- bottom collecting sump – optionally in galvanised or painted version, liquid-tight
- pull-out shelves with removable plastic sump

Cabinet for chemicals

CS-CLASSIC model CS.195.105

Body colour Anthracite grey (RAL 7016), with wing doors in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

Order No. 30495-041-30508



Cabinet for chemicals

CS-CLASSIC-G model CS.195.105.WDFW

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Pure white (RAL 9010), interior equipment with 4 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

Order No. 30496-050-30516



Cabinet for chemicals

CS-PHOENIX model CS.195.105.FD

Body colour Anthracite grey (RAL 7016), with folding doors in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

Order No. 30497-041-30508



Cabinet for chemicals


CS-CLASSIC model CS.195.105

Body colour Anthracite grey (RAL 7016), with wing doors in Gentian blue (RAL 5010), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)


Order No. 30495-047-30508




Front view



Side view



Top view



Technical data

External dimensions W x D x H
Internal dimensions W x D x H
Weight without interior equip.
Extraction air

All models

1055 x 520 x 1950 mm
996 x 493 x 1835 mm
75 kg
75 DN

Technical data

Air change **10 times**
Diff. pressure (air change **10 times**)
Maximum load
Distributed load

All models

9 m³/h
2 Pa
600 kg
421.00 kg/m²

7 door colours are available **at no extra cost**:



Order example - how to order correctly



C-LINE

Model	Information / equipment	Order No.
CS.195.105	without interior equipment	30495
CS.195.105.WDFW	without interior equipment, with wing doors incl. glass cutouts	30496
CS.195.105.FD	without interior equipment, with folding doors	30497

Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Warning yellow RAL 1004	040
Anthracite grey RAL 7016	Light grey RAL 7035	041
Anthracite grey RAL 7016	Traffic red RAL 3020	046
Anthracite grey RAL 7016	Gentian blue RAL 5010	047
Anthracite grey RAL 7016	Reseda green RAL 6011	048
Anthracite grey RAL 7016	asescos silver	049
Anthracite grey RAL 7016	Pure white RAL 9010	050

Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel galvanised	30508
3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel galvanised	30509
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated, sheet steel galvanised	30510
3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated, sheet steel galvanised	30511
3x shelf, 1x bottom collecting sump with grid (V=33.0L)	sheet steel galvanised, sheet steel painted blue	30512
3x shelf, 1x bottom collecting sump with grid (V=33.0L)	sheet steel powder-coated, sheet steel painted blue	30513
3x shelf, 1x bottom collecting sump (large) with grid (V=66.0L)	sheet steel galvanised, sheet steel painted blue	30514
3x shelf, 1x bottom collecting sump (large) with grid (V=66.0L)	sheet steel powder-coated, sheet steel painted blue	30515
4x pull-out shelf (V=33.0L)	sheet steel/polyethylene powder-coated/- grey/black	30516
5x pull-out shelf (V=33.0L)	sheet steel/polyethylene powder-coated/- grey/black	30517
6x pull-out shelf (V=33.0L)	sheet steel/polyethylene powder-coated/- grey/black	30518

Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.
shelf	sheet steel galvanised	993 x 492 x 30	100.00			30519
shelf	sheet steel powder-coated	993 x 492 x 30	100.00			30520
perforated insert	sheet steel galvanised	976 x 427 x 41	60.00			30521
bottom collecting sump	sheet steel galvanised	981 x 430 x 85*		33.00	30.00	30522
bottom collecting sump with grid	sheet steel coated blue	984 x 424 x 106*		33.00	30.00	30523
bottom collecting sump (large) with grid	sheet steel coated blue	984 x 382 x 175*		66.00	60.00	30524
pull-out shelf, with PP sump	sheet steel/polypropylene grey/black	898 x 418 x 80*	30.00	34.00	30.91	30525
box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038

Accessories	Order No.
extraction unit without exhaust air monitoring	14218

* usable storage space W x D (mm)
1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



Wolf Laboratories Limited

www.wolflabs.co.uk

Wolflabs

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

