



**Safe and approved storage of waterpollutants, non flammable liquids and toxics in working areas according to TRGS 510**

**Function / construction:**

**Cabinet body**

**robust construction and longevity:**

robust, high-quality, scratch-resistant outer body with powder-coated textured surface

**no unauthorised use:** doors

lockable with cylinder locking (integrated in the turning handle)

**easy alignment:** adjusting aids to compensate for uneven floor

**Recirculating air filter system**

- ideal for the direct installation in working areas – very low noise, only approx. 39 dB (A)
- with multi-stage wideband filter
- VDE approved monitoring electronics (extraction monitoring) incl. cool-device cable
- optical and acoustical alarm including potential-free alarm contact
- display with innovative touch panel
- manual check of the filter saturation

**Recirculating air filter storage cabinet**

**CX-CLASSIC-G model CX.229.105.WDFW**

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

**Order No. 32925-041-30508**

**Technical data**

		<b>CX.229.105.WDFW</b>
External dimensions W x D x H	mm	1055 x 520 x 2290
Internal dimensions W x D x H	mm	996 x 493 x 1835
Weight without interior equip.	kg	105
Distributed load	kg/m <sup>2</sup>	440
Frequency	Hz	50
Nominal voltage	V	230
Max. current consumption	A	0.5
Power consump. in operation	W	120
Power consump. in standby mode	W	50



**Available equipment**

- shelves – optionally in galvanised or painted version, height adjustable
- perforated insert useable as storage area
- bottom collecting sump – optionally in galvanised or painted version, liquid-tight
- pull-out shelves with removable plastic sump



**Easy operation**

Monitoring electronics with state-of-the-art 4.7" graphic display incl. touch panel, displays all essential monitoring-details



**Reusable filter system**

Multistage, high-capacity and robust broad band-pass filter, environmentally friendly due to long service lifetimes and filter recycling



**Manual test device**

Easy monitoring of the filter saturation when storing inorganic substances



**Extraction and filtration**

Filtration of solvent vapours (hydrocarbons) as well as inorganic and organic hazardous materials



**Recirculating air filter storage cabinet  
CX-CLASSIC-G model CX.229.105.WDFW**

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Pure white (RAL 9010), interior equipment with 4 x pull-out shelf (sheet steel/polyethylene powder-coated/-)


**Order No. 32925-050-30516**

**Recirculating air filter storage cabinet  
CX-CLASSIC-G model CX.229.105.WDFW**

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Light grey (RAL 7035), interior equipment with 6 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

**Order No. 32925-041-30518**

2 door colours are available **at no extra cost:**



**AVAILABLE from the first quarter of 2016**

**Order example - how to order correctly**



Model	Information / equipment	Order No.
CX.229.105.WDFW	without interior equipment, with wing doors incl. glass cutouts	32925

Body colour	Door colour	Order No.
Anthracite grey RAL 7016	Light grey RAL 7035	041
Anthracite grey RAL 7016	Pure white RAL 9010	050

Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel galvanised	30508
3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel galvanised	30509
3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated, sheet steel galvanised	30510
3x shelf, 1x bottom collecting sump (V=33.0L)	sheet steel powder-coated, sheet steel galvanised	30511
3x shelf, 1x bottom collecting sump with grid (V=33.0L)	sheet steel galvanised, sheet steel painted blue	30512
3x shelf, 1x bottom collecting sump with grid (V=33.0L)	sheet steel powder-coated, sheet steel painted blue	30513
3x shelf, 1x bottom collecting sump (large) with grid (V=66.0L)	sheet steel galvanised, sheet steel painted blue	30514
3x shelf, 1x bottom collecting sump (large) with grid (V=66.0L)	sheet steel powder-coated, sheet steel painted blue	30515
4x pull-out shelf (V=33.0L)	sheet steel/polyethylene powder-coated/- grey/black	30516
5x pull-out shelf (V=33.0L)	sheet steel/polyethylene powder-coated/- grey/black	30517
6x pull-out shelf (V=33.0L)	sheet steel/polyethylene powder-coated/- grey/black	30518

Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.
shelf	sheet steel galvanised	993 x 492 x 30	100.00			30519
shelf	sheet steel powder-coated	993 x 492 x 30	100.00			30520
perforated insert	sheet steel galvanised	976 x 427 x 41	60.00			30521
bottom collecting sump	sheet steel galvanised	981 x 430 x 85*		33.00	30.00	30522
bottom collecting sump with grid	sheet steel coated blue	984 x 424 x 106*		33.00	30.00	30523
bottom collecting sump (large) with grid	sheet steel coated blue	984 x 382 x 175*		66.00	60.00	30524
pull-out shelf, with PP sump	sheet steel/polypropylene grey/black	898 x 418 x 80*	30.00	34.00	30.91	30525
box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038

\* usable storage space W x D (mm)  
1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)

Accessories	Material	Order No.
activated carbon replacement filter	activated carbon	32952



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

