



Safe and approved storage of aggressive, non flammable hazardous materials in working areas

Function / construction:

- **robust construction and longevity:** made out of highly resistant melamin resin special plates
- **user-friendly:** access to all containers, 135° opening of the wing door, interior of the cabinet completely visible
- **no unauthorised use:** door lockable with cylinder lock and locking state indicator (red/green)
- **mobile:** castors with plinth (optionally available), easy moving of the cabinet
- **ventilation:** integrated air ducts ready for connection (DN 50) to a technical exhaust system

Available equipment

- pull-out shelves with removable plastic sump

Cabinet for acids and alkalis

SL-CLASSIC-UB model SL.080.110.UB.2T

in Light grey (RAL 7035), interior equipment with 4 x pull-out shelf (panel/PP)

Order No. 33200-001-30621



castors with plinth

stable – easy moving of the cabinet

safe fixing of the cabinet – integrated locking mechanism

flexible – height adjustment to various laboratory furniture



Cabinet for acids and alkalis

SL-CLASSIC-UB model SL.080.059.UB.TR

in Light grey (RAL 7035), interior equipment with 2 x pull-out shelf (panel/PP)

Order No. 33202-001-30625

Technical data

External dimensions W x D x H
Internal dimensions W x D x H
Weight without interior equip.
Extraction air
Air change **30 times**
Diff. pressure (air change **30 times**)
Distributed load
Maximum load

SL.080.059.UB.TR


mm 590 x 570 x 800
mm 490 x 530 x 760
kg 48
DN 50
m³/h 6
Pa 5
kg/m² 157
kg 100

SL.080.110.UB.2T


1100 x 570 x 800
all compartments: 1000 x 530 x 760
73
50
9
10
98
100




SL.080.110.UB.2T Front view




Side view




Top view




SL.080.059.UB.TR Front view



Side view



Top view



SL-LINE

Order example - how to order correctly



30623



001



30625



Model

SL.080.110.UB.2T
SL.080.059.UB.TR

Information / equipment

without interior equipment, Underbench cabinet, with 2 wing doors
without interior equipment, Underbench cabinet, right hinged, with wing door

Order No.

33200
33202



Body colour

Light grey RAL 7035

Door colour

Light grey RAL 7035

Order No.

001 4



Interior equipment packages (can only be ordered in combination with the cabinet)

Width 1100 mm
left: 1x pull-out shelf (V=11.0L)
right: 1x pull-out shelf (V=11.0L)
left: 2x pull-out shelf (V=11.0L)
right: 2x pull-out shelf (V=11.0L)

Material

panel/PP
panel/PP

Order No.

30620 2
30621 2



Interior equipment packages (can only be ordered in combination with the cabinet)

Width 590 mm
1x pull-out shelf (V=11.0L)
2x pull-out shelf (V=11.0L)

Material

panel/PP
panel/PP

Order No.

30624 2
30625 2



Interior equipment

pull-out shelf, with PP sump
castors with plinth, height = 30 mm
castors with plinth, height = 30 mm

Material

panel/PP grey
sheet steel powder-coated RAL 7035
sheet steel powder-coated RAL 7035

W x D x H (mm)

590, 1100 and 1400
1100
590



1)



2)



3)

25.00 18.00 16.36

Order No.

30770 2
30362 2
30478 2



Accessories

extraction module without exhaust air monitoring for wall mounting

Order No.

17178 2

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

