



UB-T-90 ◀ 059 cm ▶



Safety storage underbench cabinet
UB-T-90 model UB90.060.059.T

in Light grey (RAL 7035), interior equipment with 1 x perforated insert,
1 x bottom collecting sump (sheet steel powder-coated)

Order No. 30481-001-33763



Safe and approved storage of flammable hazardous substances in working areas

Function / construction:

- **robust construction and longevity:** body made of powder-coated sheet steel
- **easy handling with comfort:** door stops in any position, self-closing in the event of fire
- **user-friendly:** access to all containers, 135° opening of the wing door, interior of the cabinet completely visible
- **no unauthorised use:** door lockable with cylinder lock and locking state indicator (red/green)
- **mobile:** castors with plinth (optionally available), easy moving of the cabinet
- **ventilation:** integrated air ducts ready for connection (DN 50) to a technical exhaust system

Available equipment

- shelves with load capacity of 75 kg
- perforated insert useable as storage area
- bottom collecting sump



castors with plinth

stable – easy moving of the cabinet

safe fixing of the cabinet – integrated locking mechanism

flexible – height adjustment to various laboratory furniture

Recirculating air filter system

Avoid complex wall breakthroughs and expensive extraction piping.

Especially designed for the lower pressure conditions and volume flows of underbench cabinets.

Further information from page 395

Order No. 27144

(housing unit optional)




Technical data


		UB90.060.059.T	UB90.060.059.050.T
External dimensions W x D x H	mm	593 x 574 x 600	593 x 502 x 600
Internal dimensions W x D x H	mm	470 x 452 x 502	470 x 382 x 502
Weight without interior equip.	kg	120	110
Extraction air	DN	50	50
Air change 10 times	m ³ /h	1	1
Diff. pressure (air change 10 times)	Pa	1	1
Diff. pressure (air change 10 times EK5/AK4)	Pa	1	1
Maximum load	kg	300	300
Distributed load	kg/m ²	450.00	460.00

Usable storage heights

UB90.060.059(.050).T(R)
① mm 400*

* with perforated insert





Order example - how to order correctly



Model	Information / equipment	Order No.
UB90.060.059.T	safety storage underbench cabinet UB90.060.059.T (Light grey RAL 7035) with drawer and disposal system	30481
UB90.060.059.050.T	without interior equipment, left hinged, with wing door	30483

Body colour	Door colour	Order No.
Light grey RAL 7035	Light grey RAL 7035	001 ⁴

Depth 574 mm	Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
			old code new code
	1x perforated insert, 1x bottom collecting sump (V=12.0L)	sheet steel powder-coated RAL 7035	30485 33763 ²
	1x shelf, 1x perforated insert, 1x bottom collecting sump (V=12.0L)	sheet steel powder-coated RAL 7035	30486 33764 ²
	1x drawer (maximized internal height) (V=33.0L)	sheet steel powder-coated RAL 7035	30487 30487 ⁴
	1x perforated insert, 1x bottom collecting sump (V=11.0L)	stainless steel	30488 30488 ²
	1x shelf, 1x perforated insert, 1x bottom collecting sump (V=13.5L)	stainless steel	30489 30489 ²

Depth 502 mm	Interior equipment packages (can only be ordered in combination with the cabinet)	Material	Order No.
			old code new code
	1x perforated insert, 1x bottom collecting sump (V=10.0L)	sheet steel powder-coated RAL 7035	30490 33765 ²
	1x shelf, 1x perforated insert, 1x bottom collecting sump (V=10.0L)	sheet steel powder-coated RAL 7035	30491 33766 ²

Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.
Depth 574 mm						
shelf	sheet steel powder-coated RAL 7035	462 x 416 x 20*	75.00			30441 ²
shelf	stainless steel	462 x 416 x 20*	75.00			30442 ²
perforated insert	sheet steel powder-coated RAL 7035	420 x 439 x 60	75.00			5555 ²
perforated insert	stainless steel	425 x 439 x 60	75.00			5776 ²
bottom collecting sump	sheet steel powder-coated RAL 7035	452 x 427 x 70*		12.00	10.91	7270 ²
bottom collecting sump	stainless steel	452 x 427 x 70*		11	10	7381 ²

Interior equipment	Material	W x D x H (mm)	1)	2)	3)	Order No.
Depth 502 mm						
shelf	sheet steel powder-coated RAL 7035	461 x 344 x 20*	75.00			30368 ²
perforated insert	sheet steel powder-coated RAL 7035	442 x 345 x 70	75.00			18814 ²
bottom collecting sump	sheet steel powder-coated RAL 7035	452 x 355 x 70*		10.00	9.09	17671 ²

Interior equipment	Material	W x D x H (mm)	Order No.
All models			
pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48	30035 ²
bottle rack, for 8 bottles with Ø 90 mm	sheet steel powder-coated RAL 7035	220 x 440 x 30	7604 ²
bottle rack, for 15 bottles with Ø 70 mm	sheet steel powder-coated RAL 7035	220 x 440 x 30	7603 ²
bottle rack, for 32 bottles with Ø 40 mm	sheet steel powder-coated RAL 7035	220 x 440 x 30	7602 ²
castors with plinth, height = 30 mm	sheet steel powder-coated RAL 7035		30478 ²

Accessories	Order No.
extraction module with exhaust air monitoring for wall mounting, for underbench cabinets	24315 ²
recirculating air filter system with exhaust air monitoring, for underbench cabinets	27144 ²
adaptor, as separate housing unit for lateral mounting	26885 ²

* usable storage space W x D (mm)
 1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

