



# CAPP *Bravo*

## SINGLE CHANNEL PIPETTES

Manual pipetting has never been easier! **CAPP** has combined a lightweight ergonomic design with the latest innovations in manual liquid handling to create **CAPP *Bravo***, the most user-friendly manual pipette on the market today.

**CAPP *Bravo* Single Channel Pipette** is constructed of the highest quality autoclavable materials and features feather-light plunger and ejection forces. There is a built-in oversized indicator for easy volume monitoring and precise adjustment.

With the **CAPP *Bravo* Single Channel Pipette** all of these features are packaged in a modern, streamlined, world-renowned Danish design. Experience **CAPP**'s symphony of precise measurements and see for yourself why it is called Bravo!

### Safety:

- Fully Autoclavable
- Positive indicator lock

### Ergonomics:

- Ergonomically constructed hook for comfortable handling
- Smooth tip mounting and ejection

### Controls:

- Oversized indicator
- Reduced force feather touch plunger and ejector

### Unique Features:

- Metal shaft with double O-rings for ultra micro sizes
- Ultra lightweight construction



# Ordering Information

## CAPP *Bravo* Fixed Volume (customised volumes available on request)

| Cat. No.  | Description                            | Inaccuracy % | Imprecision % |
|-----------|--|--------------|---------------|
| B01-1F    | CAPP pipette, fixed vol. 1µL, Bravo    | 2.50         | 1.80          |
| B02-1F    | CAPP pipette, fixed vol. 2µL, Bravo    | 2.00         | 1.00          |
| B05-1F    | CAPP pipette, fixed vol. 5µL, Bravo    | 1.20         | 0.80          |
| B10-1F    | CAPP pipette, fixed vol. 10µL, Bravo   | 1.00         | 0.50          |
| B20-1F    | CAPP pipette, fixed vol. 20µL, Bravo   | 1.00         | 0.30          |
| B50-1F    | CAPP pipette, fixed vol. 50µL, Bravo   | 1.00         | 0.30          |
| B100-1F   | CAPP pipette, fixed vol. 100µL, Bravo  | 0.70         | 0.20          |
| B200-1F   | CAPP pipette, fixed vol. 200µL, Bravo  | 0.60         | 0.20          |
| B500-1F   | CAPP pipette, fixed vol. 500µL, Bravo  | 0.80         | 0.15          |
| B1000-1F  | CAPP pipette, fixed vol. 1000µL, Bravo | 0.80         | 0.20          |
| B5000-1F  | CAPP pipette, fixed vol. 5mL, Bravo    | 0.60         | 0.15          |
| B10000-1F | CAPP pipette, fixed vol. 10mL, Bravo   | 0.60         | 0.15          |

## CAPP *Bravo* Variable Volume

| Cat. No. | Description                                   | Inaccuracy % | Imprecision % | Colour |
|----------|---|--------------|---------------|--------|
| B02-1    | CAPP pipette, variable vol. 0.1-2µL, Bravo    | 12.00 / 2.00 | 6.00 / 1.00   | ■      |
| B10-1    | CAPP pipette, variable vol. 0.5-10µL, Bravo   | 2.50 / 1.00  | 1.80 / 0.50   | ○      |
| B20-1    | CAPP pipette, variable vol. 2-20µL, Bravo     | 5.00 / 1.00  | 1.50 / 0.30   | ■      |
| B50-1    | CAPP pipette, variable vol. 5-50µL, Bravo     | 2.00 / 1.00  | 1.20 / 0.30   | ■      |
| B100-1   | CAPP pipette, variable vol. 10-100µL, Bravo   | 2.50 / 0.70  | 1.00 / 0.20   | ■      |
| B200-1   | CAPP pipette, variable vol. 20-200µL, Bravo   | 2.50 / 0.60  | 0.70 / 0.20   | ■      |
| B1000-1  | CAPP pipette, variable vol. 100-1000µL, Bravo | 2.00 / 0.80  | 0.60 / 0.20   | ■      |
| B5000-1  | CAPP pipette, variable vol. 1-5mL, Bravo      | 1.20 / 0.60  | 0.30 / 0.15   | ■      |
| B10000-1 | CAPP pipette, variable vol. 1-10mL, Bravo     | 3.00 / 0.60  | 0.60 / 0.15   | ■      |

## CAPP *Bravo* Starter Packs

| Cat. No.  | Description  |
|-----------|--|
| SK-03-B01 | CAPPBravo Starter Kit including B10-1 / B100-1 / B1000-1 pipettes, C-05 stand, 3 racks of tips and 1 timer |
| SK-03-B02 | CAPPBravo Starter Kit including B10-1 / B100-1 / B1000-1 pipettes, C-05 stand, 3 racks of tips and 1 timer |
| SK-03-BB  | Three pipettes free of choice, C-05 stand, 3 racks of tips and 1 timer                                     |

## Pipette Stands

| Cat. No. | Description                        |
|----------|------------------------------------|
| C-05     | Pipette stand for up to 5 pipettes |





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

