

# ahn myLab TR-03

## Plate Rotator 30 rpm

**Usage** - ideal for processing of blood samples, separating and mixing of liquid and powdered samples

**Capacity** - 12 x 15 mL and 4 x 50 mL tubes or 22 x 15 mL tubes with the maximum operational weight capacity of 1.0 kg

**Fixed speed** - 30 rpm with an accuracy of  $\pm 2$  rpm

**Optimal mixing** - circular motion allows gentle but effective mixing

**Robust and space-saving design** - strong metal construction with rubber feet for stable operation / compact design requires minimum of bench space

**Easy maintenance and adjustment** - plates can be easily removed and installed for cleaning and replacing



## Ordering information

Description	Cat. No.
AHN myLab® TR-03 Plate Rotator 30 rpm	7-220-00-0
Plate 22 x 15 mL (included with TR-03)	7-220-02-0

Optional attachments	Cat. No.
Plate 12 x 15 mL + 4 x 50 mL	7-220-01-0

Specifications	
Run time	Continuous / infinite
Noise level	Noiseless operation
Max. loading capacity	1.0 kg
Protection class	IP 21
Dimensions (WxDxH)	346x185x145 mm
Weight	2.5 kg



Plate 22 x 15 mL



Plate 12 x 15 mL + 4 x 50 mL

### Delivery Package

1 pc. Plate Rotator 30 rpm
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card
1 pc. Plate 22 x 15 mL





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.