

Cole-Parmer®

Cole-Parmer® Handheld Homogenizer HO-200

- Ergonomic design
- Variable motor speed 5,000 to 35,000rpm
- Processing range from 0.03ml-100ml
- Can be used with robust polycarbonate disposable probes or stainless steel ones





Cole-Parmer® Handheld Homogenizer HO-200

The Cole-Parmer® HO-200 is a powerful, variable speed, homogenizer. The high-speed, high-torque motor makes the unit ideal for most homogenizing applications. During operation the rotor shaft, which is directly coupled to the motor, spins at up to 35,000rpm. The tube assembly that fits around the shaft remains static and as the rotor shaft spins, within the tube assembly, it creates a pumping action drawing the sample into the open end of the tube, or probe, and forcing it through the windows at the bottom of the tube assembly. This shearing action reduces the particle size of the sample. In addition the speed differential between the rapidly moving rotor, and the static tube assembly creates cavitation force, further dismantling the sample to additionally reduce particle size. The HO-200 is lightweight making it comfortable to use, and most samples can be processed within thirty seconds.

For longer processing the unit can be stand mounted using the included post clamp assembly. A retort stand is available separately. The HO-200 can accommodate a selection of stainless steel probes for larger and tougher samples and also our polycarbonate disposable probes. Disposable probes are suitable for hard tissues and ideal where cross contamination is of concern. If necessary disposable probes can be easily dismantled and autoclaved up to seven times. Through careful selection of probe type, samples of 0.03ml to 100ml can be processed. All Cole-Parmer® homogenizers are supplied with a tool kit for dismantling the rotor probes for easy cleaning

Technical Specification

Specification	HO-200
Motor	125 watt
Speed control	5,000 to 35,000rpm
Processing range	0.03ml to 100ml
Sound level	<72 db
Net weight, kg	0.5
Dimensions, mm (h x Ø)	160 x 55
Electrical supply	220-240V, 50Hz
IP Rating	30

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Handheld homogenizer, UK plug	04720-01	HO-200	HO-200	SHM1/UK
Handheld homogenizer, European plug	04720-02	HO-200	HO-200-EU	SHM1/EURO



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

