

# Infrared Heater IRB 1



## Description

The compact laboratory heater IRB 1 with a heating area of 60 x 60 mm is optimally suited for fast heating tasks without open flame up to max. 700°C.

### Basic equipment

Basic device with support cross

## Details

<b>Max. temperature:</b>	700°C
<b>Heating power:</b>	250 W
<b>Heating area:</b>	60 x 60 mm
<b>Max. ambient temperature:</b>	30°C

**Order no.:** 6060 000  
**Accessory:** Power control LR 1  
Support rod for stand connection

## **Technical specifications**

**Electrical supply:** 230 V, 50 / 60 Hz (115 V version on request)

**Dimensions (w x d x h):** 100 x 105 x 105 mm

**Weight:** 0,7 kg



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

