

Infrared Heater IRB 2



Description

The compact laboratory heater IRB 2 with power control and a heating area of 100 x 100 mm is optimally suited for fast heating tasks without open flame up to max. 900°C.

Basic equipment

Basic device incl. integrated power control, support cross

Details

Max. temperature:	900°C
Heating power:	800 W
Heating area:	100 x 100 mm
Power control:	stepless control between 0- 100%

Max. ambient temperature: 30°C

Order number: 6063 000

Accessory: Support rod for stand connection

Technical specifications

Electrical supply: 230 V, 50 / 60 Hz (115 V version on request)

Dimensions (w x d x h): 150 x 150 x 170 mm

Weight: 2,5 kg



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

