

Tube Roller RM 2



Description

Gentle tilting and rotating motion for uniform mixing

Basic Equipment

Basic device with 6 removable plastic rolls and spill tray

Loading Capacity

All standard tube sizes, flasks and universal vessels up to 80 mm in diameter..

Details

Motion:	rotating and tilting motion
Shaking platform:	200 x 300 mm
Max. load:	2 kg
Shaking speed:	5 - 80 rpm, in steps of 5
Angle of tilt:	3° fixed
Stroke at end of rolls:	15 mm

Dimensions of rolls: 300 x Ø 32 mm
Number of rolls: 6 pieces
Runtime: 0 - 99 h 59 min / continuous
Order No.: 6138 000
Accessory: Overhead Rack RM 2

Technical specifications

Electrical supply: 100 - 240 V, 50 / 60 Hz
Enclosure protection: IP 21
Heat emission: max. 5 W
Ambient temperature: 5°C to 45°C
Relative humidity: ~ 85 % (non-condensing)
Dimensions (w x d x h): 244 x 465 x 115 mm
Weight: 7,3 kg



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

