

Microplate Shaker TiMix 2



Description

The Microplate Shaker TiMix 2 especially for small samples in microplates, deepwell-plates or other plates in microplate format.

Basic Equipment

Basic device incl. basic platform with holders for 4 microplates

Loading Capacity

4 microplates
(deepwell-plates, or other plates in microplate format)

Details

Motion:	exact orbital motion
Shaking platform:	200 x 295 mm for 4 microplates
Max. load:	2 kg
Shaking speed:	100 - 1400 min ⁻¹ , steplessly adjustable
Stroke:	3 mm
Timer:	0 - 120 min / continuous
Order No.:	6110 000
Accessory:	Shaking platform KM Clamping strip KM

Technical specifications

Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Heat emission:	approx. 15 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	220 x 410 x 110 mm
Weight:	8,5kg



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

