


Microplate Shaker TiMix 5 control



Description

The Microplate Shaker TiMix 5 control guarantees optimal mixing of small sample volumes in microplates, deepwell-plates or other plates in microplate format. The wide choice of different rack systems allows fastening between 8 and 48 microplates. The shaker also has a digital display and programmable control.

Basic Equipment

Basic device without shaking plate, **without rack system**

With digital display and programmable control

Loading Capacity

| | |
|---|----------------|
| Standard rack system 0052 101 | 8 microplates |
| Rack system with pins 0052 096 | 8 microplates |
| Standard rack system 0052 101 + Additional tray 0052 102 | 16 microplates |
| RS high clamping pins 0052 104 | 24 microplates |
| RS metal clamping pins 0052 125 | 48 microplates |

Option:

Details

| | |
|--------------------------|--|
| Motion: | exact orbital motion |
| Shaking platform: | 400 x 300 mm for 8 microplates |
| Max. load: | 5 kg |
| Shaking speed: | 100 - 1400 rpm, in steps of 5 |
| Stroke: | 3 mm (Special stroke on request) |
| Runtime: | programmable / continuous |
| Order No.: | 6167 000 |
| Accessory: | Standard Rack System TiMix 5 Additional tray TiMix 5 Rack System with clamping pins Rack System with high clamping pins Rack System with metal clamping pins Rack System Combifix KS Universal tray KS Incubator Hood TH 15 |

Technical specifications

| | |
|--------------------------------|--|
| Electrical supply: | 230 V or 115 V, 50 / 60 H please indicate in case of order |
| Enclosure protection: | IP 21 |
| Heat emission: | approx. 25 W |
| Ambient temperature: | 5°C to 50°C |
| Relative humidity: | ~ 85 % |
| Dimensions (w x d x h): | 510 x 490 x 150 mm |
| Weight: | 20 kg |



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

