

Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines
and Vacuum Cold traps

Lablyo Extra Freezedrier

UK Manufactured

The **Lablyo Extra** offers capabilities similar to our larger H series pilot freeze driers in a small bench top system.

These include:

Control panel for programming steps using time and temperature.

Data available to download onto a memory stick.

3 temperature controlled stainless steel shelves with circulation fluid eliminates the need to pre freeze samples and allows temperature control during the freeze drying process.

4 mobile temperature probes in drying chamber for product sampling.

Condenser capacity 3kg, temperatures to -80°C

Pneumatic shelf stoppering for vial closing.



External 4 port manifold to allow for flask freeze drying.

Application:

Laboratory research and development freeze drying.

96 well plates.

Serum bottles or vials.

Bulk products for analysis.

Small R&D freeze drying.

The **Lablyo Extra** comes with clear acrylic door for observing the process.

Overall Dimensions: Height 55cm x length 70cm x width x 85cm

The **Lablyo Extra** includes a vacuum pump and connection hose that sits separately.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.