

# Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines  
and Vacuum Cold traps

## Lablyo Mini Freeze drier

### UK Manufactured

Bench top freeze drier with temperature display.  
20cm diameter x 15cm deep ice condenser  
chamber with 3kg ice capacity. 2kg ice collection  
in 24 hours.

Condenser temperature . 55°C.

Mobile temperature probe for product sampling.

### Part No's

Lablyo Mini

FIT/LYO/46/1

Vacuum pump for Lablyo

FIT/LYO/59/0

Can be used with the following freeze drying  
accessories:

### Part No's

Acrylic drying chamber

FIT/LYO/01/0

6 height shelf rack

FIT/LYO/03/0

6 heated shelf rack

FIT/LYO/19/0

4 shelf vial closing device

FIT/LYO/22/0

8 port flask manifold

FIT/LYO/13/0

48 port ampoule manifold

FIT/LYO/25/0

Spin Freezer

FIT/LYO/28/0

