

Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines
and Vacuum Cold traps

Lablyo Ultra Freezedrier

UK Manufactured

Floor standing freeze drier with new Touch Screen controls. Data logging and display of temperature and vacuum on real time graph. Downloadable via USB.
25cm diameter x 60cm deep ice condenser chamber with 10kg ice capacity. 6kg ice collection in 24 hours.
Condenser temperature -55°C .
Mobile temperature probe for product sampling.

Part No's

Lablyo Ultra FIT/LYO/56/1
Vacuum pump for Lablyo Ultra FIT/LYO/59/0

Can be used with the following freeze drying accessories:

Part No's

Acrylic drying chamber	FIT/LYO/01/0
6 height shelf rack	FIT/LYO/03/0
6 heated shelf rack	FIT/LYO/19/0
4 shelf vial closing device	FIT/LYO/22/0
8 port flask manifold	FIT/LYO/13/0
48 port ampoule manifold	FIT/LYO/25/0
Spin Freezer	FIT/LYO/28/0





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

