




VARIABLE VOLUME MECHANICAL PIPETTES



I Ultimate Performance

Built for traceability and comfort. PIPETMAN L is designed to increase daily productivity, giving you the ability to pipette longer without fatigue.

I Secure Pipetting

Reduce pipetting errors with volume-locking. Avoid accidental volume changes during pipetting cycles, ensuring your accuracy and precision.

I Optimized Comfort

Lightweight and light pipetting forces. Right-handed or left-handed, no problem. Adjust the location of the tip ejector to the most comfortable position for you.





Classic metal or plastic tip ejector

Choose between a stainless steel tip ejector for more robustness or a plastic tip ejector to avoid corrosion.



Choose your model

Choose your channel number and volume range. Each ensures consistent volume levels, allowing you to aspirate and dispense your samples with confidence.

PIPETMAN L Single Channel — Variable Volume Models				
Model	PIPETMAN DIAMOND Tips		Part Number	Volume (μL)
P2L	D10 DL10	DF10 DFL10	FA10001P FA10001M	0.2 – 2
P10L*	D10 DL10	DF10 DFL10	FA10002P FA10002M	1 – 10
P20L	D200	DF30	FA10003P FA10003M	2 – 20
P100L	D200	DF100	FA10004P FA10004M	10 – 100
P200L	D200 D300	DF200 DF300	FA10005P FA10005M	20 – 200
P1000L	D1000 D1200	DF1000 DF1200	FA10006P FA10006M	100 – 1000
P5000L	D5000		FA10007**	500 – 5000
P10mLL	D10mL		FA10008**	1000 – 10000
PIPETMAN L Multichannel				
P8x10L	D10 DL10 DF10 DFL10		FA10013	0.5 – 10
P12x10L			FA10014	
P8x20L	DL10 DFL10 DF30 D200		FA10009	2 – 20
P12x20L			FA10010	
P8x200L	D200 D300 DF200 DF300		FA10011	20 – 200
P12x200L			FA10012	
P8x300L	D200 D300 DF200 DF300		FA10015	20 – 300
P12x300L			FA10016	

For a pipette with a plastic tip ejector, the part number ends with the letter P.
For a pipette with a stainless steel tip ejector, the part number ends with the letter M.



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

