



Glenammer  
Laboratory Test Sieves

# SQ Series

## Digital Shaker



## Features

SQ digital has a digital display panel with an installed digital timer, users can set up the desired operation sieving time between 0 sec to 99 mins 59 secs

The system features a pause/start option, which enables users to operate the shaker with intervals.



### Capacity

8 X 200mm Dia. (100mm, 150mm, 3" and 8") sieves plus lid and receiver.

6 X 300mm Dia. (12" 315mm ) sieves plus lid and receiver. Extended rods available.

### Net Weight

20Kg

### Max sample weight

3Kg

### Sieve Motion

Electromagnetic

### Speed

3,000/min at 50 Hz

### Electrical Supply

230 Volt, 1 ph, 50 Hz, Input power 0.045 kw, current 0.20 amps

### Height

140mm

### Timer

Digital: 0sec - 99mins 59secs

### Net W x H x D

240mm x 140mm x 240mm

### Type

8 x 200mm dia (100mm, 150mm, 8") sieves plus lid and receiver



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

