



Glenammer
Laboratory Test Sieves

SQ Series

Analogue Shaker



Features

It has an analogue timer to adjust the operating time between 0min-60min, users can simply turn the dial to the desired time to start the operation.

SQ Analogue model support full diameter range
200 SQ Analogue is built with a stable electromagnetic motor. 200 SQ Analogue is suitable for laboratory use.
300,315, 450 SQ Analogue models are motorized shakers and are ideal for heavy-duty.



Capacity

8 X 200mm Dia. (100mm, 150mm, 3" and 8") sieves plus lid and receiver.

6 X 300mm Dia. (12" 315mm) sieves plus lid and receiver. Extended rods available.

Net Weight

17.7Kg

Max sample weight

3Kg

Sieve Motion

Electromagnetic

Speed

3,000/min at 50 Hz

Electrical Supply

230 Volt, 1 ph, 50 Hz, Input power 0.045 kw, current 0.20 amps

Height

140mm

Timer

Analogue: 0 - 60min

Net W x H x D

240mm x 140mm x 240mm

Type

8 x 200mm dia
(100mm, 150mm, 8")
sieves plus lid and receiver



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

