

## ECO MIDI K 60 CCG 4S



With an energy consumption of only 271 kWh/year (Energy test Ecodesign) the ECO MIDI 60 refrigerator ranks in energy class A.

The Ecodesign test method records the 24 hour energy consumption and multiplies it by 365 to get the annual energy consumption (AEC). The 24 hour consumption is calculated from 2 test periods of 12 hours each. The first period has sequences of door openings while the second is a stabilization period with closed door.

ECO MIDI 60 cabinets combine stylish design with high functionality and low energy consumption.

- The controls are concealed behind the top panel, protected against splashing water.
- Design principle of "no broken surfaces" makes the cabinets easy to clean.
- The distinctive aluminum profiles stretching from top to bottom on both sides of the cabinet impart strength and enhance appearance.
- The modular rails can go directly in the dishwasher.
- Easy access to stored items as the storage compartment is raised 37 cm above the floor.

ECO MIDI 60 is suitable for the 1/1 GN system. One 1/1 GN pan can be placed on each shelf.

## TECHNICAL INFORMATION

Varenummer	962600011
Category	Storage Cabinets
Type	Refrigerator
Line	Shelf 435x530 mm
Series	ECO MIDI
Equipped with	Right hand hinged reversible door, automatic door closing, 4 stainless shelves, LED lighting
Temperature range	+2/+12°C
Volume, gross	407 ltr.
Net usable volume	304 ltr.
Exterior/Interior	Stainless/Stainless
Dim. W x D x H, mm	600 x 771 x 2000
Weight, packed	114 kg
Refrigerant / GWP	R 600a (isobutane)/0.09 kg/CO2e 0.27 kg/3
Ref.capacity at -10°	275 Watt
Connection	230V, 50 Hz
Connection load	129 Watt
Energy Efficiency Class	A
Energy consumption / year (AEC)	271 kWh
Sound level	41.9 dB(A)





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

