

## SUPERIOR PLUS K 72 CCG C1 4S



With an energy consumption of only 285 kWh/year (Energy test Ecodesign) the SUPERIOR PLUS K 72 ranks in energy class A.

The Ecodesign test method records the 24 hour energy consumption and multiplies it by 365 to get the annual energy consumption (AEC). The 24 hour consumption is calculated from 2 test periods of 12 hours each. The first period has sequences of door openings while the second is a stabilization period with closed door.

- The controls are concealed behind the top panel, protected against splashing water.
- Design principle of "no broken surfaces" makes the cabinets easy to clean.
- The distinctive aluminum profiles stretching from top to bottom on both sides of the cabinet impart strength and enhance appearance.
- The modular rails can go directly in the dishwasher.
- Pedal opener with finesse - read more in the brochure.

SUPERIOR PLUS is suitable for the deep 2/1 GN system where 2 pcs. 1/1 GN pans can be placed after each other directly on the support rails.

## TECHNICAL INFORMATION

Varenummer	960720014
Category	Storage Cabinets
Type	Refrigerator
Line	2/1 GN deep
Series	SUPERIOR & ECO PLUS
Equipped with	Right hand hinged reversible door, automatic door closing, pedal door opener, 4 stainless shelves, LED lighting, castors
Temperature range	+2 / +12°C
Volume, gross	610 ltr.
Net usable volume	462 ltr.
Exterior/Interior	Stainless / Stainless
Dim. W x D x H, mm	720 x 905 x 2100
Weight, packed	136kg
Refrigerant / GWP	R 600a (isobutane)/0.07 kg/CO2e 0.22 kg /3
Ref.capacity at -10°	377 Watt
Connection	230V 50 Hz
Connection load	181 Watt
Energy Efficiency Class	A
Energy consumption / year (AEC)	285 kWh
Sound level	39.4 dB(A)





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

