


PDS No. 65580x	PRODUCT DATA SHEET			Page 1 of 1
Revision 06	96 Well Microplate, UV-Star [®] , Chimney Well			 greiner bio-one
	Greiner Item-No. 65580x			
Valid for Item-No.:	655801	655809		

1.	Description / Specification	
1.1	Description	UV-Star [®] Microplate, 96 well, clear film F-bottom (flat), alphanumeric well coding, chimney well
1.2	Dimensions	See customer drawing Foil: 135 µm (± 10 µm)
1.3	Volume	Total volume: 392 µl (mathematically calculated) Working volume: 25 - 340 µl
1.4	Material / Resin	Plate and foil: Cycloolefine, free of heavy metal
1.5	Colour	655801: plate; clear 655809: plate; black
1.6	Sterilisation	No
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	-20°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4800 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	40
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 1 December 2014	Date 2 December 2014	Date 2 December 2014	
Date 13.12.2011	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.



Wolf Laboratories Limited

www.wolflabs.co.uk

Wolflabs

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

