


Haier Biomedical

Intelligent Protection of Life Science

Wireless Monitoring Systems and Compliance Solutions



simple, smart
& easy to install



Peace of mind 24/7

Anytime, anywhere access to real time data and reports from any web enabled device.



SMS, Email or Call Alerts

Automated alerts with escalation alarms can also be set to ensure the correct people are contacted in the correct order at the correct time.



Strong, safe and robust

Utilises roaming GPRS which is easy to install and is not reliant on IT. Sensors have a signal range of 15km ensuring robust reporting and minimal infrastructure. Data is transmitted over encrypted lines to prevent hacking.



Compliant

Fully featured and documented system with training and certification. IQ/OQ/PQ available for GCLP, GMP and other areas. MKT reports for pharmaceutical storage requirements.



Flexible & Scalable



Suitable for multiple monitoring points over multiple sites – no minimum installation size. System can be expanded as client requirements change.



Automated Reports

Reports can be created and emailed to users on a weekly, monthly or yearly basis.

Wireless Monitoring Systems and Compliance Solutions



Complete Energy, Temperature and Condition Monitoring solution

A wireless real-time monitoring and control solution for energy monitoring and targeting, temperature monitoring and compliance reporting as well as condition monitoring and M2M.



Value for Money

This is an affordable and flexible product which can be retro-fitted as a standalone system or integrated with an existing BMS. Monitor a room, a building, a site or a global operation.



One System. Multiple Applications

Software modules can be tailored to a wide variety of applications with customisable dashboards and screen. For example a system purchased for one use can be extended for additional purposes without the cost of duplicating infrastructures or software packages.



No Fuss and Easy to Install

There are only three main components: wireless RF transmitters, RF gateway and Realtime Online™. No wires and no need to connect into LAN or Wi-Fi.

Wireless Monitoring Systems and Compliance Solutions

Service Contracts

Why Haier Biomedical?

- Sensor calibration traceable to UKAS 17025 standard
- System optimisation ensuring consistent reporting
- Maintain confidence in the system
- Experienced specialist engineers
- European coverage with a quick response time
- Discount on system with a large number of sensors
- Training on system operation for new and existing developments to ensure efficient system usage
- 15% discount for labour and travel costs per call out
- 10% discount for parts per call out

Haier Biomedical can also service other monitoring systems with sensor checks/calibrations.



	Standard Contract	Premium Contract
Battery strength test	✓	✓
Signal strength test	✓	✓
Recommissioning of system	✓	✓
Calibration	✓	✓
Email alarming	✓	✓
SMS alarming	✗	✓
Phone call alarming	✗	✓



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

