

## Pharmacy Refrigerators HYC-390



### Ergonomic Design



High Accuracy  
Temperature Control



Interior  
LED



Lockable  
Design



### Air Cooling, Auto Defrost

Better cabinet temperature uniformity,  
auto defrost



### Adjustable Shelves

The shelves are adjustable to any height to  
fulfill different requirements



### USB Interface

Data of temperature history can be downloaded  
via the built-in USB port and optional printer design



### Heated Double Pane Glass Door

Clear display of stored products without condensation  
at 32 °C and 85% relative humidity condition



### Self Closing Door with 90° Stay Open Feature

Energy conserving, maintains temperature  
uniformity to protect valuable samples



# Specifications

## HYC-390

### Technical Data

Cabinet Type	Upright, single glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	43

### Performance

Temp Range(℃)	2~8
---------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50/60	115/60
Power(W)	380	500
Electrical Current(A)	2.4	5.5

### Dimensions

Capacity(L/Cu.Ft)	390/13.8
Net/Gross Weight(approx)	116/139(kg) 255.7/306.4(lbs)
Interior Dimensions(W*D*H)	530*555*1380(mm) 20.9*21.9*54.3(in)
Exterior Dimensions(W*D*H)	665*710*1965(mm) 26.2*28.0*77.4(in)
Packing Dimensions(W*D*H)	735*800*2100(mm) 28.9*31.5*82.7(in)
Container Load(20'/40'/40'H)	21/42/42



### Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

### Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	7/-
USB Interface	Optional

### Others

Certificate	CE,UL
-------------	-------



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

