

Pharmacy Refrigerators HYC-610



Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Five Alarm Functions

High/low temperature, power failure, sensor error,
low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity,
auto defrost



Adjustable Shelves

The shelves are adjustable to any height to
fulfill different requirements



Heated Double Pane Glass Door

Clear display of stored products without condensation
at 32 °C and 85% relative humidity condition



Self Closing Door with 90° Stay Open Feature

Energy conserving, maintains temperature
uniformity to protect valuable samples

Specifications

HYC-610

Technical Data

Cabinet Type	Upright, single door
Climate Class	ST
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	43

Performance

Temp Range(°C)	2~8
------------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220~240/50/60	115/60
Power(W)	550	560
Electrical Current(A)	3.5	7.5

Dimensions

Capacity(L/Cu.Ft)	610/21.6
Net/Gross Weight(approx)	204/227(kg) 449.7/500.4(lbs)
Interior Dimensions(W*D*H)	680*640*1400(mm) 26.8*25.2*55.1(in)
Exterior Dimensions(W*D*H)	780*840*1945(mm) 30.7*33.1*76.6(in)
Packing Dimensions(W*D*H)	865*945*2090(mm) 34.1*37.2*82.3(in)
Container Load(20'/40'/40'H)	12/26/26



Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	6/-
USB Interface	Optional
Temp Recorder	Optional

Others

Certificate	CE
-------------	----



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

