


# Plastic Body Gel Filled ORP Electrode | HI-3230B

BNC  
connection

## Ideal for water quality control

The HI-3230B is a plastic bodied ORP/Redox probe for use with pH meters with a mV volt scale and are ideal for water quality control.

This electrode has a 1m cable and BNC connector.

- diam: 12mm
- length: 120mm
- temp: -5 to 70°C
- max pressure: 2 bar

### Recommended use

- municipal water
- quality control

### Ordering Information

#### Plastic Body Gel Filled ORP Electrode | HI-3230B

supplied with BNC connector and 1m cable

#### Electrode Storage Solutions

electrode storage solution (500 mL) HI-70300L  
electrode storage solution (230 mL) HI-70300M

### Specifications

Product Code	HI-3230B
Reference	single, Ag/AgCl
Junction	ceramic, single
Connector	BNC
Electrolyte	gel
Max Pressure	2 bar
Range	ORP: $\pm 2000$ mV; -5 to 70°C
Tip Shape	platinum pin
Temperature Sensor	No
Amplifier	No
Body Material	PEI
Recommended Use	municipal water, quality control

**BNC  
Connection**



Authorised distributor



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

