



## Great value magnetic stirrers

The HI-180 is a compact and lightweight magnetic stirrer which incorporates electronic controls to allow the user to regulate stirring speed with precision.

Equipped with the Hanna Speedsafe system, the HI-180 ensures that the safe maximum speed is never exceeded.

### Main benefits

- Neat round-edged base
- Small and compact - great for use where space is limited
- Built to last - glass-free chemical housing resists damage from accidental falls
- Safe stirring speeds guaranteed using Hanna Speedsafe technology

### Ordering Information

Magnetic stirrer | HI-180 (all models)

supplied with micro stir bar and instructions

### Choose from a great range of colours

Identify different samples or colour match your lab by choosing from nine great colour options

### Specifications

Product Code	HI-180	
Maximum Stirring Capacity	1 liter (0.26 gallons)	
Speed Range	Min.	0.01; 0.1; 1 FNU
	Max.	0.1; 1 FAU
Power Supply	110/115 VAC or 220/240 VAC, 50/60 Hz	
Cover Material	ABS plastic	
Environment	0 to 50°C (32 to 122°F) ; RH max 95%	
Dimensions / Weight	137 mm (dia) x 51 mm (h) / 640 g	





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

