

Heavy-duty Magnetic Stirrer | HI-300N



2.5 litre heavy-duty stirrer



The HI-300N is a heavy-duty magnetic stirrer with the capacity to stir up to 2.5 litres of liquid making it ideal for use in production as well as in the laboratory.

Integral Hanna Speedsafe technology ensures that the stirrer speeds can be regulated with greater precision to deliver accurate, repeatable results. The addition of a limiter promotes optimum safety by making sure that the maximum speed will not be exceeded.

Main benefits

- 2.5 litre capacity - good for use in production
- Compact, durable, efficient
- Regulate speeds with precision
- Speedsafe technology promotes optimum safety
- AISI 316 stainless steel cover - ideal for applications that create exothermic reactions


Ordering Information

Heavy duty magnetic stirrer | HI-300N

supplied with micro stir bar and instructions

Optional extra

HI-731320: extra magnetic stir bars available in packs of 10



Specifications

Product Code	HI-300N	
Maximum Stirring Capacity	2.5 litres	
Speed Range	Min.	100 rpm
	Max.	800 to 1000 rpm
Auto-Feedback	standard	
Power Supply	110/115 VAC or 230/240 VAC, 50/60 Hz	
Installation Category	II	
Cover Material	AISI 316 stainless steel	
Environment	0 to 50°C (32 to 122°F); RH max 95%	
Dimensions	180 x 180 x 70 mm	
Weight	1.4 kg	





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

