

# Cardboard Slide Folders

The cardboard slide folders by Heathrow Scientific offer a convenient solution for safely working, transporting and short-term storing of slides.

- Easily organize slides with the index on the front cover that corresponds with the number grid inside
- Recessed slide compartments and inside covers protect specimens. Slides cannot be touched when the cover is closed, reducing contamination
- Front protection flaps are attached with durable tape for repeated use. Fold flaps flat under the slide holder to save space when in use
- Thumb cut in slide compartments enables easy removal of slides with slight pressure at one end
- Compartment holder is made from a double layer of densely pressed cardboard. Cover is a single layer cardboard
- Available in two sizes that hold 20 or 30 standard (26 x 76mm and 1 x 3in) slides
- Disposable


Item No.	Color	Capacity	L x W x H in	L x W x H cm	UOM
HS9902	Black ■				
HS9913	Blue ■				
HS9914	Green ■	20 slides	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9	72/cs
HS120196	Red ■				
HS120197	Yellow ■				
HS9930	Black ■	30 slides	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1	40/cs



# Slide Folder Rack

Effortlessly locate, retrieve and keep 20-place slide folders organized

- Made from ABS
- Ideal for histology, microbiology and pathology as well as multiple user research environments
- Rack is lightweight and can be safely transported to multiple work areas
- Rack is portable and refrigerator-safe for cold sample storage
- Multiple folder racks can be stacked vertically to save space in the lab
- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached
- Additional clasps can be purchased



**Stay Organized**  
Includes 5 assorted colored clasps



HS99040A



Stackable

## Slide Folder Rack

Item No.	Color	Capacity	L x W x H in	L x W x H cm	UOM
HS99040A	Blue	5 Slide Trays	8.9 x 13.4 x 4.9	22.6 x 34 x 12.5	1/ea

*\*5 Clasps Included (Clear, Blue, Teal, Purple, Pink)*

## Slide Folder Clasps

Item No.	Color	L x W x H in	L x W x H cm	UOM
120653	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	2.2 x 1.7 x 0.6	5.6 x 4.4 x 1.6	20 ea color/bag
120654	Clear <input type="checkbox"/>			100 ea/bag

*\*Assorted colors (Clear, Gray, Blue, Green, Red)*



120653



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

