

Slide Mailers

Safe and secure transport of standard microscope slides
(76 x 26 mm and 3 x 1 in)

Cardboard Slide Mailers

Item No.	Color	Capacity	Packaging	L x W x H in	L x W x H mm
HS9903	Natural	1-Place	25	3.7 x 1.7 x 0.2	95 x 42 x 5
HS9904	Natural	2-Place	36	4.1 x 3.1 x 0.2	103 x 79 x 5
120557	Natural	4-Place	20	4.7 x 5.7 x 0.3	120 x 140 x 6

- Heavy cardboard mailers for safe transport
- Cover that opens for viewing and thumb grooves for easy slide removal
- Securely holds slides with or without cover glass
- Flat surface area is ideal for labeling or writing necessary information



Plastic Slide Mailers

Item No.	Color	Capacity	Packaging	L x W x H in	L x W x H mm
120556	Natural	1-Place	50	3.3 x 1.6 x 0.2	83 x 41 x 5
HS15983B	Natural				
HS15983A	Assorted*	2-Place	25	3.7 x 2.7 x 0.2	93 x 68 x 5

*Assorted colors (Blue, Purple, Red, Orange, Yellow)

- Polypropylene protectors ideal for secure shipping and storing of slides
- Large, easy-to-secure snaps, provide safe storage for valuable slides
- Variety of colors brings some fun to the lab and enables color coding



Plastic 5-Place Slide Mailers

Item No.	Color	Capacity	Packaging	L x W x H in	L x W x H mm
HS15982	Natural	5-Place	25	3.2 x 0.7 x 1.2	81 x 19 x 31
HS15983G	Assorted*	5-Place	25	1.7 x 0.9 x 3.5	43 x 24 x 88
HS15986	Natural				

*Assorted colors (Blue, Purple, Red, Orange, Yellow)

- Polypropylene containers ideal for secure shipping, storing or staining slides
- Polypropylene containers hold five standard slides (76 x 26 mm and 3 x 1 in)
- Inside wall features slotted channels to keep slides separated and to make removal easier
- Lid closes securely to keep slides safe



Side Opening
HS15982

End Opening
HS15986



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

