



Battery and Filter not included in warranty



Soft-Coated Triggers
with direction indicators

Ergonomically Balanced
in the hand with a loaded pipette

Universal Grip
for right- or left- handed users



Adapter
for small pipettes

RF3000®

Packed with features for comfortable and accurate pipetting

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

- Long-life li-ion battery, with quick charge—up to 20 hours use
- Toggle operation by pressing the large mode button
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes

- New USB connector ensures the appropriate charger is used with the unit
- Low-battery indicator
- Can be recharged when in use



Nose cone assembly only



Storage Stand
doubles as a charging station

Easy Identification

5 different colored nose cones help differentiate department, task or user; as well as a labeling area on the side of the unit



Magnetic & Screw Wall-Mount



3 Modes

Bright LEDs indicate 3 operation modes: high, low, & gravity



Specifications

Power	Micro-USB connector Rechargeable li-Ion battery (replaceable)
Filter	PTFE/25mm dia/0.45 µm filter
Pipettes	1–100 mL, glass and plastic
Pipette Nozzle	Polypropylene



Package Includes

- RF3000® Pipette Controller
- Charging Stand and Wall-Mountable/Magnetic Holder
- Nozzle Adapter for use with Small Pipettes
- Li-Ion Rechargeable Battery
- Five Nose Cones: Black, Yellow, Red, Grey and Blue
- 5VDC Micro USB Plug
- 6 ft (183 cm) Cable
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Item No.	Volume Range mL	L x W x H in	L x W x H cm
HS3000	1–100	9 x 8.5 x 5	22.9 x 21.6 x 12.7

Replacement parts

HS5000RFFIL	RF3000® Replacement Filters, Pack of 5		
120669	RF3000® Universal Power Adaptor 5VDC Micro USB Plug		
HS300502	RF3000® Li-Ion Rechargeable Battery, One year warranty on battery included		
120487	RF3000® Nose Cones Assorted* ■ ■ ■ ■ ■ Set of 5		
HS5000RFRFA	RF3000® Replacement Filter Adapter		
HS5000RFRP	RF3000® Replacement Pipette Adapter		

*Assorted colors (Black, Blue, Grey, Red and Yellow)



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

