



## Arctic Squares® Polycarbonate Cryostorage Boxes

Safely store vials from  $-196$  to  $+121^{\circ}\text{C}$

Our polycarbonate boxes are designed with unique features to help orient, identify and access vials with ease.

- Hinged lid with built-in stop allows easy one-handed access to samples (not on HS23325A, slip on lid)
- Forward-sloped base and high contrast indexing on the transparent lid ensure quick visual orientation
- Magnetic tipped forceps supplied with styles b and d, help pick vials and remove boxes from freezer racking


- All boxes have vent and drainage holes
- Safe for use in mechanical freezers and  $\text{LN}_2$
- Stackable
- **Autoclavable**
- US Patent No. D529622 and EP registered design No. 00551429-0001, 00551429-0002, 00551429-0003

		Part No.	Tubes	Wells	Wells dia.	L x W x H	Qty.	Color
<b>a</b>	<b>5 x 5 array boxes</b> - Ideal for small sampling - Friction fit, slide-on lid	<b>HS23325A</b>	2.0 mL	25	12.5 mm	2.99 x 2.99 x 2.0 in (7.6 x 7.6 x 5.3 cm)	8/pk	Ruby Red ■
<b>b</b>	<b>9 x 9 array standard boxes</b> - Hinged lid - Wells arranged in an easy-to-pick formation	<b>HS23381B</b>	2.0 mL	81	12.5 mm	5.23 x 5.23 x 2.0 in (13.3 x 13.3 x 5.3 cm)	4/pk	Assorted ■■■■*
<b>c</b>	<b>9 x 9 array tall boxes</b> - Use with all common sizes up to 5.0 mL - Hinged lid	<b>HS23381C</b>	5.0 mL	81	12.5 mm	5.23 x 5.23 x 3.77 in (13.3 x 13.3 x 9.6 cm)	5/pk	Amethyst Purple ■
<b>d</b>	<b>10 x 10 array boxes</b> - Maximum number of tubes in standard footprint - Hinged lid	<b>HS23300D</b>	2.0 mL	100	12.5 mm	5.23 x 5.23 x 2.0 in (13.3 x 13.3 x 5.3 cm)	4/pk	Arctic Blue ■

\*Assorted colors (Ruby Red, Amethyst, Purple, Arctic Blue, Emerald Green)

Forceps included with boxes **b**, **d**, can be used to pick vials or remove box from freezer rack.

Includes embedded magnet for easy forceps storage.





## Arctic Squares® Polycarbonate Cryostorage Boxes

Safely store vials from  $-196$  to  $+121^{\circ}\text{C}$

Our polycarbonate boxes are designed with unique features to help orient, identify and access vials with ease.

- Hinged lid with built-in stop allows easy one-handed access to samples (not on HS23325A, slip on lid)
- Forward-sloped base and high contrast indexing on the transparent lid ensure quick visual orientation
- Magnetic tipped forceps supplied with styles b and d, help pick vials and remove boxes from freezer racking
- All boxes have vent and drainage holes

- Safe for use in mechanical freezers and  $\text{LN}_2$
- Stackable
- **Autoclavable**
- US Patent No. D529622 and EP registered design No. 00551429-0001, 00551429-0002, 00551429-0003

	Part No.	Tubes	Wells	Wells Dia.	L x W x H	Qty.	Color	
<b>a</b>	<b>5 x 5 array boxes</b> - Ideal for small sampling - Friction fit, slide-on lid	<b>HS23325A</b>	2.0 mL	25	12.5 mm	2.99 x 2.99 x 2.0 in (7.6 x 7.6 x 5.3 cm)	8/pk	Ruby Red ■
<b>b</b>	<b>9 x 9 array standard boxes</b> - Hinged lid - Wells arranged in an easy-to-pick formation	<b>HS23381B</b>	2.0 mL	81	12.5 mm	5.23 x 5.23 x 2.0 in (13.3 x 13.3 x 5.3 cm)	4/pk	Assorted ■ ■ ■ ■ *
<b>c</b>	<b>9 x 9 array tall boxes</b> - Use with all common sizes up to 5.0 mL - Hinged lid	<b>HS23381C</b>	5.0 mL	81	12.5 mm	5.23 x 5.23 x 3.77 in (13.3 x 13.3 x 9.6 cm)	5/pk	Amethyst Purple ■
<b>d</b>	<b>10 x 10 array boxes</b> - Maximum number of tubes in standard footprint - Hinged lid	<b>HS23300D</b>	2.0 mL	100	12.5 mm	5.23 x 5.23 x 2.0 in (13.3 x 13.3 x 5.3 cm)	4/pk	Arctic Blue ■

\*Assorted colors (Ruby Red, Amethyst Purple, Arctic Blue, Emerald Green)



Forceps included with boxes b, d, can be used to pick vials or remove box from freezer rack

Includes embedded magnet for easy forceps storage





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

