



Reversible Tube Racks 96-Well

Easy-to-grasp tabs for sample handling and storage

96 wells for 0.5 mL tubes on one side and 96 wells for 1.5/2.0 mL tubes on the reverse side.

- Made of polypropylene
- Wells are arranged in standard 8 x 12 array
- Tabs on each end of rack snap lid into place
- 0.5 mL wells separated by molded-in grid lines, making them easier to see
- Raised rings around the 1.5 to 2.0 mL wells make tubes easier to remove
- **Autoclavable**


Part No.	Tubes	Wells	Well Dia.	L x W x H	Qty.	Color
HS2345A	1.5/2.0 mL 0.5 mL	96	11 / 7.6 mm	9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm)	5/pk	Assorted colors (Blue, Green, Pink, Yellow, Orange) *
HS2345B	1.5/2.0 mL 0.5 mL	96	11 / 7.6 mm	9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm)	5/pk	Natural
HS2345C	1.5/2.0 mL 0.5 mL	96	11 / 7.6 mm	9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm)	5/pk	Blue
HS2345D	1.5/2.0 mL 0.5 mL	96	11 / 7.6 mm	9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm)	5/pk	Green
HS2345E	1.5/2.0 mL 0.5 mL	96	11 / 7.6 mm	9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm)	5/pk	Pink
HS2345F	1.5/2.0 mL 0.5 mL	96	11 / 7.6 mm	9.6 x 4.7 x 2 in (24.6 x 12.1 x 5 cm)	5/pk	Yellow

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

Two Tier Microtube Rack 20-Well

Interlock multiple racks to increase tube storage

- Made of polypropylene
- Easily remove tubes using finger slot on front of the rack
- Racks have large molded coordinates for easy tube identification
- Bi-level rack enables easy tube label identification
- **Autoclavable**



Part No.	Tubes	Wells	Well Dia.	L x W x H	Qty.	Color
HS4284A	1.0 mL	20	10 mm	0.839 x 3.54 x 1.97 in (21.1 x 9 x 5 cm)	6/pk	Blue
	1.5 mL					
	2.0 mL					

PCR® Tube Racks 96-Well

Perfect for pre- or post-PCR® sampling

- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Tube wells are easily identified with embossed numbers and letters
- **Autoclavable**

Part No.	Tubes	Wells	Wells Dia.	L x W x H	Qty.	Color
HS23461A	0.2 mL	96	5.2 mm	5 x 3.86 x 1.2 in (13 x 9.8 x 3.3 cm)	5/pk	Assorted colors (Blue, Green, Pink, Yellow, Orange) *
HS23461B	0.2 mL	96	5.2 mm	5 x 3.86 x 1.2 in (13 x 9.8 x 3.3 cm)	5/pk	Natural
HS23461C	0.2 mL	96	5.2 mm	5 x 3.86 x 1.2 in (13 x 9.8 x 3.3 cm)	5/pk	Blue
HS23461D	0.2 mL	96	5.2 mm	5 x 3.86 x 1.2 in (13 x 9.8 x 3.3 cm)	5/pk	Green

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



PCR® is a registered trademark of Hoffman-LaRoche Inc.



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

