

4-Way Tube Rack
Link together to meet your application needs

Create as large a rack as you require for your bench top applications.

- Made of polypropylene
- Connect end to end or side to side
- Autoclavable



Part No.	Tubes	Wells	Well Dia.	L x W x H	Qty.	Color	
HS29022A	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	Natural	<input type="checkbox"/>
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				
HS29022B	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	Blue	<input checked="" type="checkbox"/>
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				
HS29022C	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	Green	<input checked="" type="checkbox"/>
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				
HS29022D	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	Pink	<input checked="" type="checkbox"/>
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				
HS29022E	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	Yellow	<input checked="" type="checkbox"/>
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				
HS29022F	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	Orange	<input checked="" type="checkbox"/>
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				
HS29022G	50 mL	4	30 mm	6.85 x 3.74 x 2.04 in (17.4 x 9.5 x 5.2 cm)	5/pk	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	*
	15 mL	12	17.5 mm				
	0.5 mL	32	8.2 mm				
	1.5/2.0 mL	32	11.6 mm				

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Cube Tube Rack
Interlocking multi-tube racks

- Made of polypropylene
- Allows you to work with different sized tubes all in one unit
- Autoclavable

Part No.	Tubes	Wells	Well Dia.	L x W x H	Qty.	Color
HS29050A	50 mL	4	30 mm	4.25 x 4.25 x 4.5 in (10.8 x 10.8 x 11.3 cm)	5/pk	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/>
	15 mL	10	17 mm			
	0.5 mL	12	11.5 mm			
	1.5/2.0 mL	16	12.7 mm			

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

