

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Reversible PCR® Rack



Removable hinged lid fits both sides

Flip this polypropylene reversible rack to the size of tubes you need. PCR side of the rack has 168 wells that hold 0.2 mL tubes. The wells accommodate 8- or 12-tube strips or individual tubes. The opposite side of the rack has 40 wells that hold 0.5 mL tubes. Both sides of the rack have 12 wells for 1.5 mL tubes. Autoclavable. Available in assorted colours: blue, green, purple, yellow, and orange.

Rack measures 206 x 131 x 54 mm.

Part Number	Description	Units
ZXZHS2344A	Reversible PCR Rack, Assorted Colours	Pck of 5



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

