


# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Website: [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

## Tube Rack 20 Well



### Interlock multiple racks to hold more tubes

Polypropylene racks lock together on each side. Easily remove tubes using finger slot on front of rack. Sturdy rack has internal ribs for stability. Racks have large molded coordinates for easy tube identification. Autoclavable. Measures 213 x 90 x 50 mm. Pack of 6.

| Part Number       | Description               | Units    |
|-------------------|---------------------------|----------|
| <b>ZXZHS4284C</b> | 20-Well Tube Racks, White | Pck of 6 |



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

