

Reax 20 - for 4 bottles (1-16 rpm)



Overhead shakers for 4 Bottles / 1 to 16 rpm

- Set and continuously adjust the variable speed on the analog control knob from 1 to 16 rpm – ideal for waste water analysis
- Easy and fast attachments of bottles
- Other vessels with 160 - 270-mm height and max. 136-mm diameter can be used with attachments

Reax 20 - for 4 bottles (1-16 rpm) - Technical Data

| | |
|--------------------------------|--|
| Overheat protection | self-resetting |
| Permissible ambient conditions | 5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity |
| Weight | 23 kg |
| Protection class DIN EN 60529 | IP 21 |
| Timer | - |
| Rotation speed range | 1 - 16 rpm |
| Power input | 280 W |
| Dimensions (w/d/h) | 490 x 520 x 465 mm |
| Max. Load | 30 kg |
| Operating Mode | - |
| Rotation speed setting | electronic control |
| Orbit / Stroke | - |
| Platform size | - |
| Accessories included | - |
| Angle | - |
| Motion | overhead |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

