

Reax 2 inc. universal adaptor



- Quickly clamps in 2 vessels of any size and accepts load capacities up to 1 kg
- Speed continuously adjustable from 20 to 100 rpm
- Universal adaptor allows for the use of various vessels of 50 - 160-mm height
- In addition, an optional adaptor is available for the use of 20 test tubes at a time

Reax 2 incl. universal adaptor - Technical Data

| | |
|--------------------------------|--|
| Overheat protection | self-resetting |
| Permissible ambient conditions | 5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity |
| Weight | 5.2 kg |
| Protection class DIN EN 60529 | IP 21 |
| Timer | - |
| Rotation speed range | 20 - 100 rpm |
| Power input | 27 W |
| Dimensions (w/d/h) | 510 x 235 x 180 mm |
| Max. Load | 1 kg |
| Operating Mode | - |
| Rotation speed setting | analog |
| Orbit / Stroke | - |
| Platform size | - |
| Accessories included | - |
| Angle | - |
| Motion | overhead |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

