



# HETTICH BENELUX SHAKERS

ENGINEERED  
FOR YOU

## Quality, Innovation and Service



is your leading supplier of laboratory equipment, which is tailored to your specific applications and requirements. Expect quality, innovation and service from Hettich Benelux, in

everything we do. You will see it in our products, A-class parts, certifications, eco-awareness, service and communication. Not only is innovation important, the maintenance of your equipment is also in good hands.

SHAKER

You can choose out of 7 different platforms, 6 different Erlenmeyer clips and further accessories which will suit your applications perfectly. It's large accessories range makes the Hettich Orbital Shaker to a real all-rounder.

All sample formats

More than 7 platforms suitable for your applications.



High temperature durability and longevity.

Solid metal housing

Angle of slope 4° / 8° as well as two platform sizes.

Application diversity

Orbital-Schüttler MO 10  
Art.-no. 982000001




- All sample formats: >7 platforms
- Solid metal housing: High temperature durability and longevity

Technical specification	
Temperature operation range	Amb. 0°C bis +65°C
Shaking frequency	50 - 400 rpm
Timer	0 - 120 min.
Maximum load	10 kg
Dimensions platform (W x H)	350 x 400 mm
Outside dimensions (W x D x H)	410 x 570 x 140 mm
Weight	13,5 kg
Electrical supply	110...230 V, 50...60 Hz

Platforms for Orbital Shaker

Double Non-slipping platform




Article no.  
300200026

Double Universal Platform with 10 rollers



Article no.  
300200028

Non-slipping platform



Article no.  
300200025

Platform for fastening Erlenmeyer clips



Article no.  
300200030

Platform for 5 separating funnels 250 ml / 5 x 500 ml / 3 x 1000 ml




Article no.  
300200037

Platform for 6 separating funnels 250 ml/ 500 ml



Article no.  
300200038

Universal Platform with 5 rollers




Article no.  
300200027

Erlenmeyer clips



Article no.	Article description	Loading capacity of platforms
300200031	Erlenmeyer clip, 50 ml	40 x 50 ml Erlenmeyer clips
300200032	Erlenmeyer clip, 100 ml	20 x 100 ml Erlenmeyer clips
300200033	Erlenmeyer clip, 250 ml	12 x 250 ml Erlenmeyer clips
300200034	Erlenmeyer clip, 500 ml	9 x 500 ml Erlenmeyer clips
300200035	Erlenmeyer clip, 1000 ml	9 x 1000 ml Erlenmeyer clips
300200036	Erlenmeyer clip 2000 ml	5 x 2000 ml Erlenmeyer clips
More accessories for Orbital Shakers		
300200029	Roller for Universal Platform	

3D Shaker 3DP 05  
Art.-no. 982200001



Technical specification	
Temperature operation range	Amb. 0°C bis +50°C
Shaking frequency	5 bis 50 rpm
Angle of slope	8°
Timer	0 - 120 min.
Maximum load	5 kg
Dimensions platform (W x H)	2 x 260 x 290 mm
Outside dimensions (W x D x H)	340 x 400 x 125 mm
Weight	8,4 kg
Electrical supply	110...230 V, 50...60 Hz

- Double capacity
- Solid metal housing: High temperature durability and longevity

More accessories for 3D Shaker	
Article no.	Article description
300200039	Rubber mat with 12 channels, 260 x 290 mm for 3D Shaker

3D Shaker 3D 05  
Art.-no. 982100001



Technical specification	
Temperature operation range	Amb. 0°C bis +50°C
Shaking frequency	5 bis 50 rpm
Angle of slope	8°
Timer	0 - 120 min.
Maximum load	5 kg
Dimensions platform (W x H)	260 x 290 mm
Outside dimensions (W x D x H)	340 x 400 x 125 mm
Weight	7,7 kg
Electrical supply	110...230 V, 50...60 Hz

- Solid metal housing
- Quiet operation! Suited for long-term work

3D Shaker 3DL 05  
Art.-no. 982300001



Technical specification	
Temperature operation range	Amb. 0°C bis +50°C
Shaking frequency	5 bis 50 rpm
Angle of slope	4°
Timer	0 - 120 min.
Maximum load	5 kg
Dimensions platform (W x H)	350 x 350 mm
Outside dimensions (W x D x H)	340 x 400 x 125 mm
Weight	8,7 kg
Electrical supply	110...230 V, 50...60 Hz

- Angle of slope 4° and extra big platform
- Solid metal housing: High temperature durability and longevity

Article no.	Article description
980000101	Country Set Europe
980000201	Country Set Switzerland
980000301	Country Set UK





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

