

IKA

designed for scientists



Hula Dancer basic

/// Data Sheet

Hula Dancer basic is a powerful orbital shaker with robust housing for high stability. The device has a shaking diameter of 4.5 mm and allows a maximum load of 100g. Its speed is fixed at 2800 rpm.



designed for scientists

Technical Data

Type of movement	orbital
Shaking stroke [mm]	4.5
Permissible shaking weight (incl. attachment) [kg]	0.1
Motor rating input [W]	9.5
Motor rating output [W]	5
Permissible ON time [%]	100
Speed fix [rpm]	2800
Timer display	none
Touch function	yes
Dimensions (W x H x D) [mm]	123 x 42 x 123
Weight [kg]	0.7
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 42
Voltage [V]	100 - 240
Frequency [Hz]	50/60
Power input standby [W]	0.7



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

