



designed for scientists



Mini MR standard

/// Data Sheet

Small magnetic stirrer without heating.

- For stirring quantities up to 1.000 ml (H₂O)
- Infinitely variable speed from 0 - 2.500 rpm
- White set-up plate suitable for observing color reactions



designed for scientists

Technical Data

Number of stirring positions	1
Stirring quantity max. per stirring position (H ₂ O) [l]	1
Motor rating output [W]	2
Direction of rotation	right
Speed control	Turning knob
Speed range [rpm]	0 - 2500
Stirring bar length [mm]	20 - 40
Set-up plate material	polyester
Set-up plate dimensions [mm]	115 x 115
Dimensions (W x H x D) [mm]	115 x 37 x 129
Weight [kg]	0.22
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 42
Voltage [V]	100 - 240
Frequency [Hz]	50/60
Power input [W]	3
DC Voltage [V=]	12
Current consumption [mA]	300



Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

