



## A500 RANGE

Solid Door Refrigerator



A500 PV



A500 PVS

### Features and Benefits

- Fan Assisted Cooling
- Auto Defrost
- Fully Automatic Operation
- Electronically Controlled
- Replaceable Door Gasket
- Digital Temperature Display
- Lock Fitted
- Reversible Door (Factory Option)
- Rollers at Rear
- Choice of White (PV) or Stainless Steel (PVS)
- Matching Freezer Page 150
- Optional Castor Base (POA)
- NB. Re-evaporation kit recommended for both high and low ambient conditions
- Adjustable Temperature

### Technical Specifications

Temperature Range	-2/+4°C
Refrigerant	R134A
Type of Defrost	Hot Gas
No of Shelves	4
Exterior Finish	White or Stainless Steel
Interior Finish	White Steel
Power Required	13 Amp Supply



Picture shows Castor Base option POA



External Dimensions	Height:	Width:	Depth:	Capacity:	Weight:
A500PV / PVS	1702mm	770mm	720mm	521ltr. (18.3cu.ft.)	94kg
Internal Dimensions					
A500PV / PVS	1494mm	650mm	555mm		

CHILLED STORAGE



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

