

## STAINLESS STEEL TROLLEYS



- Heavy Duty welded construction from 304 quality stainless steel
- Supplied in a brushed polish finish as standard
- 1.5mm thick steel worktop mounted on 40mm square tube framework
- 100mm diameter stainless steel castors with grey non-marking tyres
- Total capacity of all trolleys in the range is 150kg
- SMT3, SMTL3 & SMTL4 come with 600 x 450 x 450 single door cabinet with 1 interior shelf

Product Code	Lx W x H (mm) + caster height of x130 mm	Shelves	Weight	Format	Depth between shelves
SMT1	900 x 460 x 840	2	33KG	2 Flat Top Tiers	600 mm
SMT2	900 x 460 x 840	3	41KG	3 Flat Top Tiers	280 mm
SMT3	900 x 460 x 840	2	49KG	1 x Half Width Cabinet	600 mm
SMTL1	1200 x 550 x 840	2	44KG	2 Flat Top Tiers	600 mm
SMTL2	1200 x 550 x 840	3	52KG	3 Flat Top Tiers	280 mm
SMTL3	1200 x 550 x 840	2	60KG	1 x Half Width Cabinet	600 mm
SMTL4	1200 x 550 x 840	2	74KG	1 x Half Width Cabinet 1 x Half Width Drawer	600 mm





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

