

HAZARDOUS STORAGE
ACID & ALKALI STORAGE CABINET
88 SERIES
88A733



- High quality and reliable zintec steel cabinet for extra corrosive resistance designed to store acids safely
- This cabinet is designed for storing acids safely in its original container and is not acid proof . It is recommended that acids are stored in their original containers with the lids in the closed position.



Product Code	88A733
Hazard	Acid & Alkali
Standard Colour:	Signal White (9003)
H x W x D (mm) :	712 x 355 x 305
Aperture Size (mm):	580 x 220
Weight (kg):	14
Steel Thickness (mm):	0.9
Standard No. of Shelves:	1
Shelf Size (mm):	345 x 230 x 25
Shelf Capacity (ltrs):	2
Additional Shelf Code:	88ES33A
Number of Doors:	1— Reinforced, rebated
Lock:	3 way locking chrome handle
Accessories Included:	Adjustable spillage tray
Shelf Capacity UDL (kg):	75
Certification:	DSEAR 2015 HS (G) 51 1990 Factory Inspectorate Certificate No 1 Parts 3&4
Hazard label:	White diamond, red border, black test tubes, black hand
Lead Time:	15 days + delivery

KINDLY NOTE A SUMP TRAY IS NOT PROVIDED WITH THIS SIZE CABINET



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

