

HAZARDOUS STORAGE FLAMMABLE STORAGE CABINET

88 SERIES

88F794



- High quality and reliable storage cabinet designed to maintain integrity for 30 minutes during fire
- Creates a safe and secure storage solution for flammable and hazardous materials when stored according to the material's own storage guidelines



Product Code	88F794
Hazard	Flammable
Standard Colour:	Golden Yellow (BS08E51)
H x W x D (mm):	712 x 915 x 459
Aperture Size (mm):	580 x 779
Weight (kg):	35
Steel Thickness (mm):	0.9
Standard No. of Shelves:	1
Shelf Size (mm):	900 x 385 x 40
Shelf Capacity (ltrs) :	13
Additional Shelf Code:	88ES94F
Optional Floorstand:	88FFS3 - 450 x 915 x 459
Number of Doors:	2 — reinforced, rebated
Lock:	3 way locking chrome handle
Accessories Included:	Adjustable spillage tray Liquid tight sump tray
Shelf Capacity UDL (kg):	75
Certification:	DSEAR 2015 HS (G) 51 1990 Factory Inspectorate Certificate No 1 Parts 3&4
Hazard label:	White diamond, red border, black flame
Sump Tray Capacity (ltrs):	31
Sump Tray Size (mm):	770 x 410 x 100
Lead Time:	15 days + delivery



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

