

HAZARDOUS STORAGE
STOCK FLAMMABLE STORAGE CABINET
88 SERIES
88F794STOCK



- High quality and reliable storage cabinet designed to maintain integrity for 30 minutes during fire
- Creates a safe and secure storage solution for flammable and hazardous materials when stored according to the material's



| Product Code | 88F794STOCK |
|--|---|
| Hazard | Flammable |
| Standard Colour: | Golden Yellow (BS08E51) |
| H x W x D (mm): | 712 x 915 x 459 |
| Aperture Size (mm): | 562 x 830 |
| Weight (kg): | 38 |
| Steel Thickness (mm): | 0.9 |
| Standard No. of Shelves: | 1 |
| Please note all shelves are perforated | |
| Shelf Size (mm): | 910 x 410 x 30 |
| Additional Shelf Code: | 88ES94FSTOCK |
| Optional Floorstand: | 88FFS3 - 450 x 915 x 459 |
| Number of Doors: | 2 — reinforced, rebated |
| Lock: | 2 way locking chrome handle |
| Accessories Included: | Adjustable shelf Liquid tight sump tray |
| Shelf Capacity UDL (kg): | 75 |
| Certification: | DSEAR 2015 HS (G) 51 1990 Factory Inspectorate Certificate No 1 Parts 3&4 |
| Hazard label: | Red diamond, black border, black flame |
| Sump Tray Capacity (ltrs): | 20 |
| Sump Tray Size (mm): | 820 X 410 x 80 |
| Lead Time: | Immediate despatch |



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

