

**HAZARDOUS STORAGE**  
**ACID & ALKALI STORAGE CABINET**  
**88 SERIES**  
**88A644**



- High quality and reliable zintec steel cabinet for extra corrosive resistance designed to store acids safely
- This cabinet is designed for storing acid safely in its original container and is not acid proof. It is recommended that acids are stored in their original containers with the lids in the closed position.



Product Code	88A644
Hazard	Acid & Alkali
Standard Colour:	Signal White (9003)
H x W x D (mm):	610 x 459 x 459
Aperture Size (mm):	475 x 320
Weight (kg):	18
Steel Thickness (mm):	0.9
Standard No. of Shelves:	1
Shelf Size (mm):	450 x 385 x 40
Shelf Capacity (ltrs) :	6.5
Additional Shelf Code:	88ES44A
Optional Floorstand:	88AFS2 450 x 459 x 459
Number of Doors:	1— reinforced, rebated
Lock:	3 way locking chrome handle
Accessories Included:	Adjustable spillage tray Liquid tight sump tray
Shelf Capacity UDL (kg):	75
Certification:	DSEAR 2015 HS (G) 51 1990 Factory Inspectorate Certificate No 1 Parts 3&4
Hazard label:	White diamond, red border, black test tubes, black hand
Sump Tray Capacity (ltrs):	12.5
Sump Tray Size (mm):	310 x 410 x 100
Lead Time:	15 days + delivery



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

